

Recombinant Human Luteinizing Hormone

货号: **AYB27230**

Recombinant Human Luteinizing Hormone

产品信息

生物活性	Recombinant Human Luteinizing Hormone is a recombinant form of Luteinizing hormone (human) (HY-P2 293). Luteinizing hormone (human), a heterodimeric glycoprotein hormone produced by the pituitary gland (LH), plays key roles in human reproduction.
CAS	
中文名称	
分子量	
体外研究	
体内研究	
形式	Solid
运输条件	Room temperature in continental US; may vary elsewhere.
保存条件	Please store the product under the recommended conditions in the Certificate of Analysis.
溶解性	
纯度	