


# POPSO

货号: **AYB27370**

## 产品信息

生物活性	POPSO is a <b>zwitterionic buffer</b> , increases osmolality and shows marked inhibition of anion uniport. POPSO inhibits chloride uniport with an <b>IC<sub>50</sub></b> value of 24 mM. POPSO enhances copper uptake and toxicity in alga, impairs mitochondrial inner membrane. The working pH range of POPSO sesquisodium salt is 7.2-8.5.
CAS	68189-43-5
中文名称	
分子量	362.42
体外研究	
体内研究	
形式	Solid
运输条件	Room temperature in continental US; may vary elsewhere.
保存条件	
溶解性	
纯度	