

Triosephosphate isomerase

货号: **AYB27421**

Triosephosphate isomerase

产品信息

| | |
|------|---|
| 生物活性 | Triosephosphate isomerase (TPI) is a glycolytic enzyme. Triosephosphate isomerase fast interconverts dihydroxyacetone phosphate and D: -glyceraldehyde-3-phosphate, its catalytic site is at the dimer interface. |
| CAS | 9023-78-3 |
| 中文名称 | |
| 分子量 | |
| 体外研究 | |
| 体内研究 | |
| 形式 | |
| 运输条件 | Room temperature in continental US; may vary elsewhere. |
| 保存条件 | Please store the product under the recommended conditions in the Certificate of Analysis. |
| 溶解性 | |
| 纯度 | |