


DOTAP Transfection Reagent

货号: **AYB27429**

产品信息

生物活性	DOTAP Transfection Reagent is a cationic derivative of trimethylammonium attached to two 18-carbon fatty acid tails, each with a single alkene group. 18:1 TAP is a cationic liposome-forming compound useful for the transfection of DNA, RNA, and other negatively charged molecules into eukaryotic cells. Reagent grade, for research use only.
CAS	144189-73-1
中文名称	(2,3-二油酰基-丙基)-三甲胺(MS盐)DOTAP(methyl sulfate salt)
分子量	774.19
体外研究	
体内研究	
形式	
运输条件	Room temperature in continental US; may vary elsewhere.
保存条件	Please store the product under the recommended conditions in the Certificate of Analysis.
溶解性	
纯度	