


1,2-Dioleoyl-sn-glycerol

货号: **AYB27448**

产品信息

生物活性	ALN29882 is a glycerolipid located on the plasma membrane. It consists of two fatty acid chains covalently linked to a single glycerol molecule by means of an ester bond. 18:1 DG has been used as a source of diacylglycerol in the diacylglycerol O-acyltransferase 1 (DGAT1) assay. It is also used as a substrate in the DGAT-1 enzyme assay to evaluate compounds as potential inhibitors of DGAT-1. Suitable for lipoprotein overlay screening assays with the recombinant protein His-AtROP6.
CAS	24529-88-2
中文名称	(S)-1,2-二油精
分子量	620.99
体外研究	
体内研究	
形式	
运输条件	Room temperature in continental US; may vary elsewhere.
保存条件	Please store the product under the recommended conditions in the Certificate of Analysis.
溶解性	
纯度	