


Lumigen APS-5

货号: **AYB27481**



产品信息

生物活性	Lumigen APS-5 is a substrate of alkaline phosphatase (ALP) . Lumigen APS-5 can be used to assess the activity of alkaline phosphatase (ALP) .
CAS	193884-53-6
中文名称	APS-5 (4-氯苯巯基)(10-甲基-9,10-二氢吡啶亚甲基)磷酸二钠盐
分子量	491.84
体外研究	<p>Guidelines (Following is our recommended protocol. This protocol only provides a guideline, and should be modified according to your specific needs).</p> <p>碱性磷酸酶活性测定:</p> <ol style="list-style-type: none"> 1. 将测试样品 (10 mg) 与 Lumigen APS-5 (100 μL) 混合, 并在最高设置下旋转 10 秒。对照组加入 10 μL 蒸馏水。 2. 混合样品在 75$^{\circ}$C 水浴中热处理 1 分钟以灭活非细菌 ALP。热处理后, 样品放置在 96 孔板中立即测量。用光度计测定 ALP 的活性。 <p>The accuracy of these methods have not been independently confirmed. They are for reference only.</p>
体内研究	
形式	
运输条件	Room temperature in continental US; may vary elsewhere.
保存条件	Please store the product under the recommended conditions in the Certificate of Analysis.
溶解性	
纯度	