


TAPS

货号: **AYB27502**

产品信息

生物活性	TAPS is a biological buffer, remain lysozyme native structure intact and prevents thermal denaturation against high temperatures. TAPS exhibits pK _a value of 8.1, while the half-maximum values of connexin channel activity is 8.5 (pH).
CAS	29915-38-6
中文名称	
分子量	243.28
体外研究	
体内研究	
形式	Solid
运输条件	Room temperature in continental US; may vary elsewhere.
保存条件	
溶解性	
纯度	