


# TES sodium

货号: **AYB27504**

## 产品信息

生物活性	TES sodium is a <b>buffering agent</b> ( $pK_a=7.550$ at $25^\circ\text{C}$ ). TES sodium is one of the Good's buffers, the buffer capacity ranging pH 6.8-8.2.
CAS	70331-82-7
中文名称	
分子量	251.23
体外研究	
体内研究	
形式	
运输条件	Room temperature in continental US; may vary elsewhere.
保存条件	Please store the product under the recommended conditions in the Certificate of Analysis.
溶解性	
纯度	