


4-Methylumbelliferyl- β -D-glucuronide

货号: **AYB27512**



产品信息

生物活性	4-Methylumbelliferyl- β -D-glucuronide is a fluorogenic substrat (λ_{ex} =362 nm, λ_{em} =445 nm). 4-Methylumbelliferyl- β -D-glucuronide has potential applications in detecting the activity of β -glucuronidase and the number of <i>Escherichia coli</i> .
CAS	6160-80-1
中文名称	4-甲基-2-氧代-2H-1-苯并吡喃-7-基-beta-D-葡糖苷酸
分子量	352.29
体外研究	4-Methylumbelliferyl- β -D-glucuronide (156 μ M-5 mM; 30 min) 在 37°C 下释放 4-甲基伞形酮, 用于酶促反应。 4-Methylumbelliferyl- β -D-glucuronide (50 μ g/mL; 24 h) 在 TSA 培养基中用于检测大肠杆菌数。 The accuracy of these methods have not been independently confirmed. They are for reference only.
体内研究	
形式	
运输条件	Room temperature in continental US; may vary elsewhere.
保存条件	Please store the product under the recommended conditions in the Certificate of Analysis.
溶解性	
纯度	