

Adenylate

货号: **AYB27532**

Adenylate

产品信息

生物活性	Adenylate (Adk), i.e., myokinase, is a phosphotransferase enzyme, is often used in biochemical studies. Adenylate catalyzes the interconversion of adenosine phosphates. Adenylate monitors phosphate nucleotide levels inside the cell, it plays an important role in cellular energy homeostasis.
CAS	9013-02-9
中文名称	
分子量	
体外研究	
体内研究	
形式	
运输条件	Room temperature in continental US; may vary elsewhere.
保存条件	Please store the product under the recommended conditions in the Certificate of Analysis.
溶解性	
纯度	