


# Guanosine 5'-triphosphate-5'-adenosine

货号: **AYB27547**



## 产品信息

生物活性	Guanosine 5'-triphosphate-5'-adenosine, the 5' <b>cap analog</b> , is a fluorescent substrate analog.
<b>CAS</b>	10527-47-6
中文名称	
分子量	772.41
体外研究	<p>Guanosine 5'-triphosphate-5'-adenosine (GpppA) labeled with pyrene at the 3'-O position of adenosine acts as an artificial substrate.</p> <p>Fluorescently labeled Guanosine 5'-triphosphate-5'-adenosine (GpppA) and GpppG analogs as potential substrates that represent a reasonable compromise between the structural complexity and requirements of the enzyme.</p> <p><b>The accuracy of these methods have not been independently confirmed. They are for reference only.</b></p>
体内研究	
形式	
运输条件	Room temperature in continental US; may vary elsewhere.
保存条件	Please store the product under the recommended conditions in the Certificate of Analysis.
溶解性	
纯度	