


# Hexanoic anhydride

货号: **AYB27679**

## 产品信息

生物活性	Hexanoic anhydride is used as a reactant for the synthesis of cremomannolipin A. For the green synthesis of acyclovir dipivoxil (acyclovir prodrug), for the preparation of hexanoyl-modified chitosan nanoparticles, chitosan-based polymer surfactants by N-acylation of chitosan .
CAS	2051-49-2
中文名称	正己酸酐
分子量	214.30
体外研究	
体内研究	
形式	
运输条件	Room temperature in continental US; may vary elsewhere.
保存条件	Please store the product under the recommended conditions in the Certificate of Analysis.
溶解性	
纯度	