


BES sodium

货号: **AYB27724**

产品信息

生物活性	BES sodium is a secondary standard biochemical buffer with effective pH range of 6.4 to 7.8 (pKa: 7.15 at 20 °C). BES sodium is used in the diagnostic testing area.
CAS	66992-27-6
中文名称	
分子量	235.23
体外研究	
体内研究	
形式	Solid
运输条件	Room temperature in continental US; may vary elsewhere.
保存条件	4°C, sealed storage, away from moisture
溶解性	
纯度	