


# Cyanine 3 Bihexanoic Acid Dye, Potassium Salt

货号: **AYB27782**

## 产品信息

生物活性	Cyanine 3 Bihexanoic Acid Dye, Potassium Salt is a biochemical reagent that can be used as a biological material or organic compound for life science related research.
CAS	762260-71-9
中文名称	
分子量	755.96
体外研究	
体内研究	
形式	
运输条件	Room temperature in continental US; may vary elsewhere.
保存条件	Please store the product under the recommended conditions in the Certificate of Analysis.
溶解性	
纯度	