


MES monohydrate

货号: **AYB27820**

产品信息

生物活性	MES (2-Morpholinoethanesulphonic acid) monohydrate is a buffering agent in biology and biochemistry . MES monohydrate is one of the Good's buffers, the buffer capacity ranging pH 5.5-7.0. MES monohydrate is broadly used to regulate pH value for plants culture medium, reagent solution, and physiological experiments.
CAS	145224-94-8
中文名称	MES(一水合物)
分子量	213.25
体外研究	
体内研究	
形式	Solid
运输条件	Room temperature in continental US; may vary elsewhere.
保存条件	4°C, sealed storage, away from moisture
溶解性	
纯度	