

eGFP mRNA-LNP

货号: **AYB27866**

eGFP mRNA-LNP

产品信息

生物活性	eGFP mRNA-LNP is an enhanced green fluorescent protein. After eGFP mRNA enters the cells, green fluorescent protein will be expressed. And lipid nanoparticles (LNP) encapsulated with eGFP mRNA are suitable for detection of RNA delivery, translation efficiency, cell viability, etc.
CAS	
中文名称	
分子量	
体外研究	
体内研究	
形式	
运输条件	Room temperature in continental US; may vary elsewhere.
保存条件	Please store the product under the recommended conditions in the Certificate of Analysis.
溶解性	
纯度	