


Phenylethyl resorcinol

货号: **AYB27870**

产品信息

| | |
|------|--|
| 生物活性 | Phenylethyl resorcinol is a biochemical reagent that can be used as a biological material or organic compound for life science related research. |
| CAS | 85-27-8 |
| 中文名称 | 4-(α -甲基苯甲基)间苯二酚 |
| 分子量 | 214.26 |
| 体外研究 | |
| 体内研究 | |
| 形式 | Solid |
| 运输条件 | Room temperature in continental US; may vary elsewhere. |
| 保存条件 | 4°C, stored under nitrogen |
| 溶解性 | |
| 纯度 | |