

Angiotensin-converting enzyme

货号: **AYB27958**

Angiotensin-converting enzyme

产品信息

生物活性	Angiotensin-converting enzyme (ACE) is a dicarboxypeptidase, it converts inactive Angiotensin I (Ang I) to active Ang II and degrades active bradykinin (BK). Angiotensin-converting enzyme is a potent vasoconstrictor, is often used in biochemical studies.
CAS	9015-82-1
中文名称	
分子量	
体外研究	
体内研究	
形式	
运输条件	Room temperature in continental US; may vary elsewhere.
保存条件	Please store the product under the recommended conditions in the Certificate of Analysis.
溶解性	
纯度	