

D-Mannose isomerase

货号: **AYB27996**

D-Mannose isomerase

产品信息

| | |
|------|--|
| 生物活性 | D-Mannose isomerase is a carbohydrate isomerase catalyzing the reversible isomerization of d-mannose to d-fructose. D-Mannose isomerase belongs to the N-acetylglucosamine 2-epimerase (AGE) superfamily along with AGE, cellobiose 2-epimerase (CE), and aldose-ketose isomerase (AKI). |
| CAS | 9031-25-8 |
| 中文名称 | |
| 分子量 | |
| 体外研究 | |
| 体内研究 | |
| 形式 | |
| 运输条件 | Room temperature in continental US; may vary elsewhere. |
| 保存条件 | Please store the product under the recommended conditions in the Certificate of Analysis. |
| 溶解性 | |
| 纯度 | |