

CD77 synthase

货号: **AYB28045**

CD77 synthase

产品信息

生物活性	CD77 synthase (A4GALT) is a glycosphingolipid-specific glycosyltransferase. CD77 synthase transfers a galactose to the alpha-1,4 position of lactosylceramide to form globotriaosylceramide. CD77 synthase can be used for the synthesis of P1 blood group antigens.
CAS	52725-57-2
中文名称	
分子量	
体外研究	
体内研究	
形式	
运输条件	Room temperature in continental US; may vary elsewhere.
保存条件	Please store the product under the recommended conditions in the Certificate of Analysis.
溶解性	
纯度	