


# MES potassium

货号: **AYB28068**

## 产品信息

生物活性	MES (2-Morpholinoethanesulphonic acid) potassium is a <b>buffering agent</b> in biology and biochemistry. MES potassium is one of the Good's buffers, the buffer capacity ranging pH 5.5-7.0. MES potassium is broadly used to regulate pH value for <b>plants</b> culture medium, reagent solution, and physiological experiments.
CAS	39946-25-3
中文名称	吗啉乙磺酸钾盐
分子量	233.33
体外研究	
体内研究	
形式	
运输条件	Room temperature in continental US; may vary elsewhere.
保存条件	Please store the product under the recommended conditions in the Certificate of Analysis.
溶解性	
纯度	