


ZnAF-2

货号: **AYB28109**



产品信息

生物活性	ZnAF-2, 6-iso- is a cell-impermeable compound that acts as a high-affinity Zn ²⁺ -specific fluorescent probe (K _d = 2.7 nM). It exhibits low basal fluorescence with an approximately 51-fold increase in fluorescence intensity upon stoichiometric (1:1) binding to Zn ²⁺ . Little affinity for Ca ²⁺ , Mg ²⁺ , Na ⁺ or K ⁺ .
CAS	321859-11-4
中文名称	
分子量	572.61
体外研究	
体内研究	
形式	
运输条件	Room temperature in continental US; may vary elsewhere.
保存条件	Please store the product under the recommended conditions in the Certificate of Analysis.
溶解性	
纯度	