

# Carboxypeptidase C

货号: AYB28126

## Carboxypeptidase C

### 产品信息

生物活性	<b>Carboxypeptidase C</b> is a <b>carboxypeptidase</b> , is often used in biochemical studies. <b>Carboxypeptidase C</b> removes COOH-terminal lysine, arginine, and proline, as well as all other neutral, aliphatic, aromatic, and the acidic protein <b>amino acids</b> of a peptide chain.
CAS	9046-67-7
中文名称	
分子量	
体外研究	Carboxypeptidase C (羧肽酶 C) 在 pH 4.5-6.0 时活性最佳。 <b>The accuracy of these methods have not been independently confirmed. They are for reference only.</b>
体内研究	
形式	
运输条件	Room temperature in continental US; may vary elsewhere.
保存条件	Please store the product under the recommended conditions in the Certificate of Analysis.
溶解性	
纯度	