


Triarachidonin

货号: **AYB28175**

产品信息

生物活性	Triarachidonin is a polyunsaturated triacylglycerol with arachidonoyl (20:4) side chains at the C-1, C-2, and C-3 positions, useful in lipid research . N/A
CAS	23314-57-0
中文名称	三酰肌醇
分子量	951.45
体外研究	
体内研究	
形式	
运输条件	Room temperature in continental US; may vary elsewhere.
保存条件	Please store the product under the recommended conditions in the Certificate of Analysis.
溶解性	
纯度	