

# Inorganic pyrophosphatase

货号: **AYB28266**

## Inorganic pyrophosphatase

### 产品信息

生物活性	Inorganic pyrophosphatase (PPase) converts pyrophosphate (PPi) to phosphate. Inorganic pyrophosphatase is an essential component of in vitro transcription reactions for RNA preparation, is often used in biochemical studies. Inorganic pyrophosphatase is critical for driving cellular processes such as nucleic acid and protein synthesis.
CAS	9024-82-2
中文名称	
分子量	
体外研究	
体内研究	
形式	
运输条件	Room temperature in continental US; may vary elsewhere.
保存条件	Please store the product under the recommended conditions in the Certificate of Analysis.
溶解性	
纯度	