

Phosphoglucomutase

货号: **AYB28281**

Phosphoglucomutase

产品信息

生物活性	Phosphoglucomutase is often used in biochemical studies. Phosphoglucomutase is an enzyme that can transfer the phosphate group on the α -D-glucose monomer forward from the 1-position to the 6-position or reversely transfer from the 6-position to the 1-position, and promote the glucose-1-phosphate and glucose-6-phosphate Transform each other. Phosphoglucomutase is a key enzyme in glycolysis and gluconeogenesis, and plays an important role in the metabolism of proteins, lipids and nucleic acids.
CAS	9001-81-4
中文名称	
分子量	
体外研究	
体内研究	
形式	
运输条件	Room temperature in continental US; may vary elsewhere.
保存条件	Please store the product under the recommended conditions in the Certificate of Analysis.
溶解性	
纯度	