


# 5-Phenyl-1H-pyrrol-3-yl tosyl-L-alaninate

货号: **AYB28323**

## 产品信息

生物活性	5-Phenyl-1H-pyrrol-3-yl tosyl-L-alaninate is a biochemical reagent that can be used as a biological material or organic compound for life science related research.
CAS	99740-00-8
中文名称	吡咯酯3-(N-TOSYL-L-ALANINYLAZY)-5-PHENYLPYRROLE
分子量	384.45
体外研究	
体内研究	
形式	
运输条件	Room temperature in continental US; may vary elsewhere.
保存条件	Please store the product under the recommended conditions in the Certificate of Analysis.
溶解性	
纯度	