


# 18:0 EPC chloride

货号: **AYB28336**

## 产品信息

|      |  |
|------|--|
| 生物活性 | 18:0 EPC chloride is a synthetic cationic phospholipid. 18:0 EPC chloride (at the critical synergistic concentrations of 2.34-2.93 $\mu\text{M}$ ) significantly improves the inactivation effect of eugenol against <i>Escherichia coli</i> . |
| CAS  | 328268-13-9  |
| 中文名称 |  |
| 分子量  | 854.66   |
| 体外研究 |  |
| 体内研究 |  |
| 形式   |  |
| 运输条件 | Room temperature in continental US; may vary elsewhere.  |
| 保存条件 | Please store the product under the recommended conditions in the Certificate of Analysis.  |
| 溶解性  |  |
| 纯度   |  |