


cis-13,16,19-Docosatrienoic acid methyl ester

货号: **AYB28357**

产品信息

生物活性	cis-13,16,19-Docosatrienoic acid methyl ester is an ester form of docosatrienoic acid, a rare omega-3 fatty acid not found in the normal phospholipid polyunsaturated fatty acid pool. Easily detectable, it inhibits the binding of [3H]-LTB4 to porcine neutrophils at a concentration of 5 μ M. Methyl docosatrienoate is more fat-soluble than the free acid and may be more desirable than the free acid in some formulations.
CAS	108698-01-7
中文名称	(13Z,16Z,19Z)-13,16,19-二十二碳三烯酸甲酯
分子量	348.56
体外研究	
体内研究	
形式	
运输条件	Room temperature in continental US; may vary elsewhere.
保存条件	Please store the product under the recommended conditions in the Certificate of Analysis.
溶解性	
纯度	