

Vaccinia virus capping enzyme

货号: **AYB28359**

Vaccinia virus capping enzyme

产品信息

生物活性	Vaccinia virus capping enzyme is a transcription initiation factor. Vaccinia virus capping enzyme is a heterodimer of D1 (844 aa) and D12 (287 aa) polypeptides that executes all three steps in m7GpppRNA synthesis. Vaccinia virus capping enzyme has been used widely as a reagent for capping and cap-labeling RNAs <i>in vitro</i> .
CAS	
中文名称	
分子量	
体外研究	Vaccinia virus capping enzyme 参与 5' 帽结构的形成和早期转录的终止。 The accuracy of these methods have not been independently confirmed. They are for reference only.
体内研究	
形式	
运输条件	Room temperature in continental US; may vary elsewhere.
保存条件	Please store the product under the recommended conditions in the Certificate of Analysis.
溶解性	
纯度	