


Poly(ethylene glycol) tetrahydrofurfuryl ether

货号: **AYB28365**

产品信息

| | |
|------|--|
| 生物活性 | Poly(ethylene glycol) tetrahydrofurfuryl ether is liquid glycolen can be easily transformed into a gel system with excellent elasticity, so it can be used as a medium for dissolving water-insoluble drugs. |
| CAS | 31692-85-0 |
| 中文名称 | 二甘醇单四氢糠基醚和三甘醇单四氢糠基醚的混合物 |
| 分子量 | |
| 体外研究 | |
| 体内研究 | |
| 形式 | |
| 运输条件 | Room temperature in continental US; may vary elsewhere. |
| 保存条件 | Please store the product under the recommended conditions in the Certificate of Analysis. |
| 溶解性 | |
| 纯度 | |