


TAPSO

货号: **AYB28374**

产品信息

生物活性	TAPSO is a common biological buffer, a pH stabilization reagent, with effective pH ranging from 7-8, and pKa values ranging from 7.5-9.0. TAPSO contains Tris groups and exhibits quite reactive activity with zwitterionic glycine peptides .
CAS	68399-81-5
中文名称	
分子量	259.28
体外研究	
体内研究	
形式	Solid
运输条件	Room temperature in continental US; may vary elsewhere.
保存条件	4°C, protect from light
溶解性	
纯度	