

— ABLYBIO, Help Your Research



Interferon alpha/beta Receptor 1 (YD12293) Rabbit mAb

货号: **AYD11023**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB FC IP
推荐浓度	
理论分子量	64kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SP2/0,HeLa,A-549
细胞定位	Cell membrane, Late endosome, Lysosome
纯化	亲和纯化


抗原信息

抗原信息	
------	--

靶点信息

研究背景	The protein encoded by this gene is a type I membrane protein that forms one of the two chains of a receptor for interferons alpha and beta. Binding and activation of the receptor stimulates Janus protein kinases, which in turn phosphorylate several proteins, including STAT1 and STAT2. The encoded protein also functions as an antiviral factor.
基因ID	3454
基因名	IFNAR1
Swiss	P17181 (https://www.uniprot.org/uniprotkb/P17181/entry)
别名	Interferon alpha/beta Receptor 1 (YD12293), Interferon alpha/beta Receptor 1 (YD12293) Rabbit mAb, IFNAR1, Cytokine receptor class-II member 1, Cytokine receptor family 2 member 1, Type I interferon receptor 1, IFNAR

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)