

— ABLYBIO, Help Your Research



# n-Myc (YD11489) Rabbit mAb

货号: **AYD11043**

## 产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IP
推荐浓度	
理论分子量	50kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA, 50% glycerol, pH7.3.
偶联物	Unconjugated
阳性对照	HeLa
细胞定位	Nucleus
纯化	亲和纯化

## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	This gene is a member of the MYC family and encodes a protein with a basic helix-loop-helix (bHLH) domain. This protein is located in the nucleus and must dimerize with another bHLH protein in order to bind DNA. Amplification of this gene is associated with a variety of tumors, most notably neuroblastomas. Multiple alternatively spliced transcript variants encoding different isoforms have been found for this gene.
------	---

基因ID	4613
基因名	MYCN
Swiss	P04198 ( <a href="https://www.uniprot.org/uniprotkb/P04198/entry">https://www.uniprot.org/uniprotkb/P04198/entry</a> )
别名	n-Myc (YD11489),n-Myc (YD11489) Rabbit mAb,MYCN,Class E basic helix-loop-helix protein 37,BHLHE37,NMYC

### 产品验证

### 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)