

— ABLYBIO, Help Your Research



beta III Tubulin (YD14072) Rabbit mAb

货号: **AYD11074**

产品信息

| | |
|-------|---|
| 反应 | Human,Mouse,Rat |
| 宿主 | Rabbit |
| 克隆性 | Monoclonal |
| 预测反应 | |
| 应用 | WB IHC-P ICC/IF FC IP |
| 推荐浓度 | |
| 理论分子量 | 50kDa |
| 实测分子量 | |
| 形式 | Liquid |
| 保存条件 | Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3. |
| 偶联物 | Unconjugated |
| 阳性对照 | HeLa,293T,MCF7,mouse brain,mouse testis,rat brain,rat testis |
| 细胞定位 | Cytoplasm, cytoskeleton, Cell projection, growth cone, lamellipodium, filopodium |
| 纯化 | 亲和纯化 |


抗原信息

| | |
|------|--|
| 抗原信息 | |
|------|--|

靶点信息

| | |
|-------|---|
| 研究背景 | This gene encodes a class III member of the beta tubulin protein family. Beta tubulins are one of two core protein families (alpha and beta tubulins) that heterodimerize and assemble to form microtubules. This protein is primarily expressed in neurons and may be involved in neurogenesis and axon guidance and maintenance. Mutations in this gene are the cause of congenital fibrosis of the extraocular muscles type 3. Alternate splicing results in multiple transcript variants. A pseudogene of this gene is found on chromosome 6. |
| 基因ID | 10381 |
| 基因名 | TUBB3 |
| Swiss | Q13509 (https://www.uniprot.org/uniprotkb/Q13509/entry) |
| 别名 | beta III Tubulin (YD14072),beta III Tubulin (YD14072) Rabbit mAb,TUBB3,Tubulin beta-4 chain,Tubulin beta-III,TUBB4 |

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)