

— ABLYBIO, Help Your Research



# Cullin 3 (YD12872) Rabbit mAb

货号: **AYD11131**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB FC
推荐浓度	
理论分子量	89kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-251MG,HeLa,MCF7,Mouse brain,Mouse spleen,Mouse testis,Rat brain
细胞定位	Nucleus, Golgi apparatus, Cell projection, cilium, flagellum, Cytoplasm, cytoskeleton, spindle, microtubule organizing center, centrosome, spindle pole
纯化	亲和纯化




## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	This gene encodes a member of the cullin protein family. The encoded protein plays a critical role in the polyubiquitination and subsequent degradation of specific protein substrates as the core component and scaffold protein of an E3 ubiquitin ligase complex. Complexes including the encoded protein may also play a role in late endosome maturation. Mutations in this gene are a cause of type 2E pseudohypoaldosteronism. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.
基因ID	8452
基因名	CUL3
Swiss	Q13618 ( <a href="https://www.uniprot.org/uniprotkb/Q13618/entry">https://www.uniprot.org/uniprotkb/Q13618/entry</a> )
别名	Cullin 3 (YD12872),Cullin 3 (YD12872) Rabbit mAb,CUL3,KIAA0617

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)