

— ABLYBIO, Help Your Research



# Nucleolin (YD32633) Rabbit mAb

货号: **AYD11144**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	
理论分子量	77kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	MCF7,A-549,THP-1
细胞定位	Nucleus, nucleolus, Cytoplasm
纯化	亲和纯化

## 抗原信息


抗原信息	
------	--

## 靶点信息

研究背景	Nucleolin (NCL), a eukaryotic nucleolar phosphoprotein, is involved in the synthesis and maturation of ribosomes. It is located mainly in dense fibrillar regions of the nucleolus. Human NCL gene consists of 14 exons with 13 introns and spans approximately 11kb. The intron 11 of the NCL gene encodes a small nucleolar RNA, termed U20.
------	------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

基因ID	4691
基因名	NCL
Swiss	P19338 ( <a href="https://www.uniprot.org/uniprotkb/P19338/entry">https://www.uniprot.org/uniprotkb/P19338/entry</a> )
别名	Nucleolin (YD32633),Nucleolin (YD32633) Rabbit mAb,NCL,Protein C23

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)