

— ABLYBIO, Help Your Research



# TBK1 (YD12313) Rabbit mAb

货号: **AYD11183**

## 产品信息

反应	Human,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P IP
推荐浓度	
理论分子量	84kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,HeLa,Mouse lung,Rat thymus
细胞定位	Cytoplasm
纯化	亲和纯化


## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	The NF-kappa-B (NFKB) complex of proteins is inhibited by I-kappa-B (IKB) proteins, which inactivate NFKB by trapping it in the cytoplasm. Phosphorylation of serine residues on the IKB proteins by IKB kinases marks them for destruction via the ubiquitination pathway, thereby allowing activation and nuclear translocation of the NFKB complex. The protein encoded by this gene is similar to IKB kinases and can mediate NFKB activation in response to certain growth factors.
基因ID	29110
基因名	TBK1
Swiss	Q9UHD2 ( <a href="https://www.uniprot.org/uniprotkb/Q9UHD2/entry">https://www.uniprot.org/uniprotkb/Q9UHD2/entry</a> )
别名	TBK1 (YD12313),TBK1 (YD12313) Rabbit mAb,TBK1,NF-kappa-B-activating kinase,T2K,TANK-binding kinase 1,NAK

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)