

— ABLYBIO, Help Your Research



# TRPV5 (YD19572) Rabbit mAb

货号: **AYD11189**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P
推荐浓度	
理论分子量	83kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,293F,Mouse kidney,Mouse heart,Rat heart,Rat kidney
细胞定位	Apical cell membrane
纯化	亲和纯化

## 抗原信息


抗原信息	
------	--

## 靶点信息

研究背景	This gene is a member of the transient receptor family and the TrpV subfamily. The calcium-selective channel encoded by this gene has 6 transmembrane-spanning domains, multiple potential phosphorylation sites, an N-linked glycosylation site, and 5 ANK repeats. This protein forms homotetramers or heterotetramers and is activated by a low internal calcium level.
------	--

基因ID	56302
基因名	TRPV5
Swiss	Q9NQA5 ( <a href="https://www.uniprot.org/uniprotkb/Q9NQA5/entry">https://www.uniprot.org/uniprotkb/Q9NQA5/entry</a> )
别名	TRPV5 (YD19572),TRPV5 (YD19572) Rabbit mAb,TRPV5,Calcium transport protein 2,Epithelial calcium channel 1,Osm-9-like TRP channel 3,ECAC1

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)