

Ferritin Heavy Chain?(YD13104) Rabbit mAb

货号: **AYD11194**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB ICC/IF FC
推荐浓度	
理论分子量	21kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse heart
细胞定位	Cytoplasm, Lysosome, Cytoplasmic vesicle, autophagosome
纯化	

抗原信息





抗原信息	
------	--

靶点信息

研究背景	This gene encodes the heavy subunit of ferritin, the major intracellular iron storage protein in prokaryotes and eukaryotes. It is composed of 24 subunits of the heavy and light ferritin chains. Variation in ferritin subunit composition may affect the rates of iron uptake and release in different tissues. A major function of ferritin is the storage of iron in a soluble and nontoxic state. Defects in ferritin proteins are associated with several neurodegenerative diseases. This gene has multiple pseudogenes. Several alternatively spliced transcript variants have been observed, but their biological validity has not been determined.
基因ID	2495
基因名	FTH1

Swiss	P02794
别名	Ferritin Heavy Chain?(YD13104)

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn