

— ABLYBIO, Help Your Research



Follistatin (YD19527) Rabbit mAb

货号: **AYD11232**

产品信息

反应	Human,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	
理论分子量	38kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,Rat pancreas
细胞定位	Secreted, Nucleus, nucleolus
纯化	亲和纯化


抗原信息

抗原信息	
------	--

靶点信息

研究背景	Follistatin is a single-chain gonadal protein that specifically inhibits follicle-stimulating hormone release. The single FST gene encodes two isoforms, FST317 and FST344 containing 317 and 344 amino acids respectively, resulting from alternative splicing of the precursor mRNA. In a study in which 37 candidate genes were tested for linkage and association with polycystic ovary syndrome (PCOS) or hyperandrogenemia in 150 families, evidence was found for linkage between PCOS and follistatin.
基因ID	10468
基因名	FST
Swiss	P19883 (https://www.uniprot.org/uniprotkb/P19883/entry)
别名	Follistatin (YD19527),Follistatin (YD19527) Rabbit mAb,FST,Activin-binding protein

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)