

— ABLYBIO, Help Your Research



SCARB1 (YD19552) Rabbit mAb

货号: **AYD11235**

产品信息

| | |
|-------|---|
| 反应 | Human,Mouse,Rat |
| 宿主 | Rabbit |
| 克隆性 | Monoclonal |
| 预测反应 | |
| 应用 | WB IHC-P |
| 推荐浓度 | |
| 理论分子量 | 61kDa |
| 实测分子量 | |
| 形式 | Liquid |
| 保存条件 | Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3. |
| 偶联物 | Unconjugated |
| 阳性对照 | 293T,Rat ovary |
| 细胞定位 | Cell membrane, Membrane, caveola |
| 纯化 | 亲和纯化 |

抗原信息



| | |
|------|--|
| 抗原信息 | |
|------|--|

靶点信息

| | |
|------|---|
| 研究背景 | The protein encoded by this gene is a plasma membrane receptor for high density lipoprotein cholesterol (HDL). The encoded protein mediates cholesterol transfer to and from HDL. In addition, this protein is a receptor for hepatitis C virus glycoprotein E2. Two transcript variants encoding different isoforms have been found for this gene. |
|------|---|

| | |
|-------|---|
| 基因ID | 949 |
| 基因名 | SCARB1 |
| Swiss | Q8WTV0 (https://www.uniprot.org/uniprotkb/Q8WTV0/entry) |
| 别名 | SCARB1 (YD19552),SCARB1 (YD19552) Rabbit mAb,SCARB1,CD36 and LIMPII analogous 1,CD36 antigen-like 1,Collagen type I receptor,thrombospondin receptor-like 1,SR-BI,CD36L1,CLA1 |

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)