

— ABLYBIO, Help Your Research



# CXCR5 (YD34355) Rabbit mAb

货号: **AYD11300**

## 产品信息

反应	Human, Mouse, Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC ICC/IF FC
推荐浓度	
理论分子量	42kDa/42kDa/42kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA, 50% glycerol, pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Cell membrane
纯化	亲和纯化


## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	Cytokine receptor that binds to B-lymphocyte chemoattractant (BLC). Involved in B-cell migration into B-cell follicles of spleen and Peyer patches but not into those of mesenteric or peripheral lymph nodes. May have a regulatory function in Burkitt lymphoma (BL) lymphomagenesis and/or B-cell differentiation Cytokine receptor that binds to B-lymphocyte chemoattractant (BLC). Involved in B-cell migration into B-cell follicles of spleen and Peyer patches but not into those of mesenteric or peripheral lymph nodes Cytokine receptor that binds to B-lymphocyte chemoattractant (BLC). Involved in B-cell migration into B-cell follicles of spleen and Peyer patches but not into those of mesenteric or peripheral lymph nodes (By similarity)
基因ID	643
基因名	CXCR5, Cxcr5
Swiss	P32302 ( <a href="https://www.uniprot.org/uniprotkb/P32302/entry">https://www.uniprot.org/uniprotkb/P32302/entry</a> ), Q04683 ( <a href="https://www.uniprot.org/uniprotkb/Q04683/entry">https://www.uniprot.org/uniprotkb/Q04683/entry</a> ), P34997 ( <a href="https://www.uniprot.org/uniprotkb/P34997/entry">https://www.uniprot.org/uniprotkb/P34997/entry</a> )
别名	CXCR5 (YD34355),CXCR5 (YD34355) Rabbit mAb,CXCR5,Burkitt lymphoma receptor 1,Monocyte-derived receptor 15,Burkitt lymphoma receptor 1 homolog,Neurolymphatic receptor,BLR1,MDR15,Gpcr6

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)