

— ABLYBIO, Help Your Research



DAB2 (YD21087) Rabbit mAb

货号: **AYD11306**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P IP
推荐浓度	
理论分子量	82kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa
细胞定位	Cytoplasm, Cytoplasmic vesicle, clathrin-coated vesicle membrane, Membrane, clathrin-coated pit
纯化	亲和纯化



抗原信息

抗原信息	
------	--

靶点信息

研究背景	This gene encodes a mitogen-responsive phosphoprotein. It is expressed in normal ovarian epithelial cells, but is down-regulated or absent from ovarian carcinoma cell lines, suggesting its role as a tumor suppressor. This protein binds to the SH3 domains of GRB2, an adaptor protein that couples tyrosine kinase receptors to SOS (a guanine nucleotide exchange factor for Ras), via its C-terminal proline-rich sequences, and may thus modulate growth factor/Ras pathways by competing with SOS for binding to GRB2. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.
基因ID	1601
基因名	DAB2
Swiss	P98082 (https://www.uniprot.org/uniprotkb/P98082/entry)
别名	DAB2 (YD21087),DAB2 (YD21087) Rabbit mAb,DAB2,Adaptor molecule disabled-2,Differentially expressed in ovarian carcinoma 2,Differentially-expressed protein 2,DOC2

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)