

# DAB2 (YD21087) Rabbit mAb

货号: **AYD11306**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P IP
推荐浓度	
理论分子量	82kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa
细胞定位	Cytoplasm, Cytoplasmic vesicle, clathrin-coated vesicle membrane, Membrane, clathrin-coated pit
纯化	

## 抗原信息



抗原信息	
------	--

## 靶点信息

研究背景	This gene encodes a mitogen-responsive phosphoprotein. It is expressed in normal ovarian epithelial cells , but is down-regulated or absent from ovarian carcinoma cell lines, suggesting its role as a tumor suppressor. This protein binds to the SH3 domains of GRB2, an adaptor protein that couples tyrosine kinase receptors to SOS (a guanine nucleotide exchange factor for Ras), via its C-terminal proline-rich sequences, and may thus modulate growth factor/Ras pathways by competing with SOS for binding to GRB2. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.
基因ID	1601
基因名	DAB2

Swiss	P98082
别名	DAB2 (YD21087)

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)