

— ABLYBIO, Help Your Research



RAGE (YD19309) Rabbit mAb

货号: **AYD11336**

产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	
理论分子量	43kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Rat lung
细胞定位	Cell membrane, Cell projection, phagocytic cup, Early endosome, Nucleus, Secreted
纯化	亲和纯化


抗原信息

抗原信息	
------	--

靶点信息

研究背景	The advanced glycosylation end product (AGE) receptor encoded by this gene is a member of the immunoglobulin superfamily of cell surface receptors. It is a multiligand receptor, and besides AGE, interacts with other molecules implicated in homeostasis, development, and inflammation, and certain diseases, such as diabetes and Alzheimer's disease. Many alternatively spliced transcript variants encoding different isoforms, as well as non-protein-coding variants, have been described for this gene (PMID:18089847).
基因ID	177
基因名	AGER
Swiss	Q15109 (https://www.uniprot.org/uniprotkb/Q15109/entry)
别名	RAGE (YD19309),RAGE (YD19309) Rabbit mAb,AGER,Receptor for advanced glycosylation end products,RAGE

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)