



# ENO1/ENO2/ENO3 (YD12372) Rabbit mAb

货号: **AYD11354**

产品信息

|       |                                                                                               |
|-------|-----------------------------------------------------------------------------------------------|
| 反应    | Human,Mouse,Rat                                                                               |
| 宿主    | Rabbit                                                                                        |
| 克隆性   | Monoclonal                                                                                    |
| 预测反应  |                                                                                               |
| 应用    | WB ICC/IF FC                                                                                  |
| 推荐浓度  |                                                                                               |
| 理论分子量 | 47kDa                                                                                         |
| 实测分子量 |                                                                                               |
| 形式    | Liquid                                                                                        |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated                                                                                  |
| 阳性对照  | HepG2,SW620,A-549,HeLa,Mouse heart,Mouse skeletal muscle,Rat heart                            |
| 细胞定位  | Cytoplasm                                                                                     |
| 纯化    | 亲和纯化                                                                                          |


抗原信息

|      |  |
|------|--|
| 抗原信息 |  |
|------|--|

靶点信息

|       |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
|-------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景  | This gene encodes one of the three enolase isoenzymes found in mammals. This isoenzyme is found in skeletal muscle cells in the adult where it may play a role in muscle development and regeneration. A switch from alpha enolase to beta enolase occurs in muscle tissue during development in rodents. Mutations in this gene have been associated with glycogen storage disease. Alternatively spliced transcript variants encoding different isoforms have been described. |
| 基因ID  | 2027                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| 基因名   | ENO3                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| Swiss | P13929 ( <a href="https://www.uniprot.org/uniprotkb/P13929/entry">https://www.uniprot.org/uniprotkb/P13929/entry</a> )                                                                                                                                                                                                                                                                                                                                                          |
| 别名    | ENO1/ENO2/ENO3 (YD12372), ENO1/ENO2/ENO3 (YD12372) Rabbit mAb, ENO3, 2-phospho-D-glycerate hydro-lyase, Enolase 3, Muscle-specific enolase, Skeletal muscle enolase                                                                                                                                                                                                                                                                                                             |

产品验证



实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)