

IL-26 (YD19685) Rabbit mAb

货号: **AYD11376**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P ICC/IF IP
推荐浓度	
理论分子量	20kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse heart
细胞定位	Secreted
纯化	

抗原信息



抗原信息	
------	--

靶点信息

研究背景	This gene was identified by its overexpression specifically in herpesvirus samimiri-transformed T cells. The encoded protein is a member of the IL10 family of cytokines. It is a secreted protein and may function as a homodimer. This protein is thought to contribute to the transformed phenotype of T cells after infection by herpesvirus samimiri. [provided by RefSeq, Jul 2008]
基因ID	55801
基因名	IL26
Swiss	Q9NPH9

别名	IL-26 (YD19685)
----	-----------------

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn