

— ABLYBIO, Help Your Research



# Phospho-PAK4 (Ser474)/PAK5 (Ser602)/PAK6 (Ser560) (YD32661) Rabbit mAb

货号: AYD11380

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB FC
推荐浓度	
理论分子量	64kDa/75kDa/81kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Cytoplasm, Nucleus, Mitochondrion
纯化	亲和纯化


## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	<p>Serine/threonine-protein kinase that plays a role in a variety of different signaling pathways including cytoskeleton regulation, cell adhesion turnover, cell migration, growth, proliferation or cell survival (PubMed:26598620). Activation by various effectors including growth factor receptors or active CDC42 and RAC1 results in a conformational change and a subsequent autophosphorylation on several serine and/or threonine residues. Phosphorylates and inactivates the protein phosphatase SSH1, leading to increased inhibitor y phosphorylation of the actin binding/depolymerizing factor cofilin. Decreased cofilin activity may lead to stabilization of actin filaments. Phosphorylates LIMK1, a kinase that also inhibits the activity of cofilin. Phosphorylates integrin beta5/ITGB5 and thus regulates cell motility. Phosphorylates ARHGEF2 and activates the downstream target RHOA that plays a role in the regulation of assembly of focal adhesions and actin stress fibers. Stimulates cell survival by phosphorylating the BCL2 antagonist of cell death BAD. Alternatively, inhibits apoptosis by preventing caspase-8 binding to death domain receptors in a kinase independent manner. Plays a role in cell-cycle progression by controlling levels of the cell-cycle regulatory protein CDKN1A and by phosphorylating RAN. Promotes kinase-independent stabilization of RHOA, thereby contributing to focal adhesion disassembly during cell migration (PubMed:26598620)</p> <p>Serine/threonine protein kinase that plays a role in the regulation of gene transcription. The kinase activity is induced by various effectors including AR or MAP2K6/MAPKK6. Phosphorylates the DNA-binding domain of androgen receptor/AR and thereby inhibits AR-mediated transcription. Also inhibits ESR1-mediated transcription. May play a role in cytoskeleton regulation by interacting with IQGAP1. May protect cells from apoptosis through phosphorylation of BAD</p> <p>Serine/threonine protein kinase that plays a role in a variety of different signaling pathways including cytoskeleton regulation, cell migration, proliferation or cell survival. Activation by various effectors including growth factor receptors or active CDC42 and RAC1 results in a conformational change and a subsequent autophosphorylation on several serine and/or threonine residues. Phosphorylates the proto-oncogene RAF1 and stimulates its kinase activity. Promotes cell survival by phosphorylating the BCL2 antagonist of cell death BAD. Phosphorylates CTNND1, probably to regulate cytoskeletal organization and cell morphology. Keeps microtubules stable through MARK2 inhibition and destabilizes the F-actin network leading to the disappearance of stress fibers and focal adhesions</p>
基因ID	10298, 106821730, 57144
基因名	PAK4, PAK6, PAK5
Swiss	O96013 ( <a href="https://www.uniprot.org/uniprotkb/O96013/entry">https://www.uniprot.org/uniprotkb/O96013/entry</a> ), Q9NQU5 ( <a href="https://www.uniprot.org/uniprotkb/Q9NQU5/entry">https://www.uniprot.org/uniprotkb/Q9NQU5/entry</a> ), Q9P286 ( <a href="https://www.uniprot.org/uniprotkb/Q9P286/entry">https://www.uniprot.org/uniprotkb/Q9P286/entry</a> )
别名	Phospho-PAK4 (Ser474)/PAK5 (Ser602)/PAK6 (Ser560) (YD32661),Phospho-PAK4 (Ser474)/PAK5 (Ser602)/PAK6 (Ser560) (YD32661) Rabbit mAb,PAK4,PAK6,PAK5,p21-activated kinase 4,PAK-5,p21-activated kinase 6,p21-activated kinase 5,p21-activated kinase 7,KIAA1142,KIAA1264,PAK7

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/>[www.ablybio.cn](http://www.ablybio.cn))