

— ABLYBIO, Help Your Research



# Poliovirus Receptor (YD33088) Rabbit mAb

货号: **AYD11384**

## 产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC-P
推荐浓度	
理论分子量	45kDa
实测分子量	
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-549,Daudi(negative)
细胞定位	Cell membrane, Secreted
纯化	亲和纯化



## 抗原信息

抗原信息	
------	--

## 靶点信息

研究背景	The protein encoded by this gene is a transmembrane glycoprotein belonging to the immunoglobulin superfamily. The external domain mediates cell attachment to the extracellular matrix molecule vitronectin, while its intracellular domain interacts with the dynein light chain Tctex-1/DYNLT1. The gene is specific to the primate lineage, and serves as a cellular receptor for poliovirus in the first step of poliovirus replication. Multiple transcript variants encoding different isoforms have been found for this gene.
基因ID	5817
基因名	PVR
Swiss	P15151 ( <a href="https://www.uniprot.org/uniprotkb/P15151/entry">https://www.uniprot.org/uniprotkb/P15151/entry</a> )
别名	Poliovirus Receptor (YD33088), Poliovirus Receptor (YD33088) Rabbit mAb, PVR, Nectin-like protein 5, PVS

## 产品验证



## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)